

For Research Use Only

ACIII Recombinant antibody, PBS Only (Capture)

Catalog Number: 84492-4-PBS



Basic Information

Catalog Number:

84492-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004036

GeneID (NCBI):

109

UNIPROT ID:

O60266

Full Name:

adenylate cyclase 3

Calculated MW:

129 kDa

Purification Method:

Protein A purification

CloneNo.:

241906E9

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01346-1, ADCY3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84492-4-PBS. Detection antibody: 84492-2-PBS. Standard: Ag25859. Range: 1.563-100 ng/mL.