

For Research Use Only

TTL12 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84493-2-PBS



Basic Information

Catalog Number: 84493-2-PBS	GenBank Accession Number: BC001070	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23170	CloneNo.: 241904D7
Source: Rabbit	UNIPROT ID: Q14166	
Isotype: IgG	Full Name: tubulin tyrosine ligase-like family, member 12	
Immunogen Catalog Number: AG13804	Calculated MW: 644 aa, 74 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01347-1, TTLL12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84493-2-PBS. Detection antibody: 84493-3-PBS. Standard: Ag13804. Range: 0.625-80 ng/mL