

For Research Use Only

TTL12 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84493-3-PBS



Basic Information

Catalog Number:

84493-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13804

GenBank Accession Number:

BC001070

GeneID (NCBI):

23170

UNIPROT ID:

Q14166

Full Name:

tubulin tyrosine ligase-like family, member 12

Calculated MW:

644 aa, 74 kDa

Purification Method:

Protein A purification

CloneNo.:

241904E3

Applications

Tested Applications:

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01347-1, TTLL12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84493-2-PBS. Detection antibody: 84493-3-PBS. Standard: Ag13804. Range: 0.625-80 ng/mL