

For Research Use Only

LMO7 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84495-3-PBS



Basic Information

Catalog Number: 84495-3-PBS	GenBank Accession Number: NM_005358	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4008	CloneNo.: 241781E11
Source: Rabbit	UNIPROT ID: Q8WWI1	
Isotype: IgG	Full Name: LIM domain 7	
Immunogen Catalog Number: AG30016	Calculated MW: 193 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01336-1, LMO7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84495-1-PBS. Detection antibody: 84495-3-PBS. Standard: Ag30016. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01336-3, LMO7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84495-2-PBS. Detection antibody: 84495-3-PBS. Standard: Ag30016. Range: 0.625-40 ng/mL