

For Research Use Only

IL-2 R gamma/CD132 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84501-4-PBS



Basic Information

| | | |
|---------------------------------------|--|---|
| Catalog Number: 84501-4-PBS | GenBank Accession Number: NM_000206.3 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 3561 | CloneNo.: 242020G10 |
| Source: Rabbit | UNIPROT ID: P31785-1 | |
| Isotype: IgG | Full Name: interleukin 2 receptor, gamma (severe combined immunodeficiency) | |
| | Calculated MW: NM_000206.3 | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01367-1, IL-2RG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84501-2-PBS. Detection antibody: 84501-4-PBS. Standard: Eg1705. Range: 0.313-40 ng/mL.