

For Research Use Only

Mouse RBP4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84503-2-PBS



Basic Information

Catalog Number:

84503-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_011255.3

GeneID (NCBI):

19662

UNIPROT ID:

Q00724

Full Name:

retinol binding protein 4, plasma

Calculated MW:

23kDa

Purification Method:

Protein A purification

CloneNo.:

242040D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01365-2, MOUSE RBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84503-1-PBS. Detection antibody: 84503-2-PBS. Standard: Eg1140. Range: 0.156-10 ng/mL