

For Research Use Only

CD3E Recombinant antibody, PBS Only (Capture)

Catalog Number: 84505-3-PBS



Basic Information

Catalog Number: 84505-3-PBS	GenBank Accession Number: BC049847	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 916	CloneNo.: 242012G3
Source: Rabbit	ENSEMBL Gene ID: ENSG00000198851	
Isotype: IgG	UNIPROT ID: P07766	
	Full Name: CD3e molecule, epsilon (CD3-TCR complex)	
	Calculated MW: 207 aa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01376-2, CD3E Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84505-3-PBS. Detection antibody: 84505-2-PBS. Standard: Eg2144. Range: 0.313-40 ng/mL