

For Research Use Only

CD37 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84509-3-PBS



Basic Information

| | | |
|---------------------------------------|---|---|
| Catalog Number: 84509-3-PBS | GenBank Accession Number: NM_001774.3 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 951 | CloneNo.: 242011B12 |
| Source: Rabbit | UNIPROT ID: P11049-1 | |
| Isotype: IgG | Full Name: CD37 molecule | |
| | Calculated MW: 32kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01362-1, CD37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84509-3-PBS. Detection antibody: 84509-2-PBS. Standard: Eg1764. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01362-2, CD37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84509-3-PBS. Detection antibody: 84509-1-PBS. Standard: Eg1764. Range: 0.313-40 ng/mL