

For Research Use Only

TMEM41B Recombinant antibody, PBS Only (Capture)

Catalog Number: 84510-3-PBS



Basic Information

Catalog Number: 84510-3-PBS	GenBank Accession Number: BC035034	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 440026	CloneNo.: 242017E3
Source: Rabbit	UNIPROT ID: Q5BJD5	
Isotype: IgG	Full Name: transmembrane protein 41B	
Immunogen Catalog Number: AG30810		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01380-1, TMEM41B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84510-3-PBS. Detection antibody: 84510-1-PBS. Standard: Ag30810. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01380-2, TMEM41B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84510-3-PBS. Detection antibody: 84510-2-PBS. Standard: Ag30810. Range: 0.625-80 ng/mL