

For Research Use Only

TMEM41B Recombinant antibody

Catalog Number: 84510-4-RR



Basic Information

Catalog Number: 84510-4-RR	GenBank Accession Number: BC035034	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 440026	CloneNo.: 242017C11
Source: Rabbit	UNIPROT ID: Q5BJD5	Recommended Dilutions: WB 1:500-1:2000
Isotype: IgG	Full Name: transmembrane protein 41B	
Immunogen Catalog Number: AG30810	Calculated MW: 32 kDa	
	Observed MW: 25-28 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HEK-293 cells, A549 cells, A2780 cells
Species Specificity: human	

Background Information

Transmembrane protein 41B (TMEM41B), also known as Stasimon, is a major target of survival motor neuron (SMN)-dependent U12 splicing and has been identified as a protein required for motor circuit function (PMID: 23063131). TMEM41B localizes to the endoplasmic reticulum. It is involved in autophagosome biogenesis and lipid mobilization (PMID: 30126924). TMEM41B has been shown to be essential for mouse embryonic development (PMID: 30352685). Additionally, genome-scale CRISPR knockout screens reveal that TMEM41B is required for infection of flavivirus and coronavirus, including SARS-CoV-2 (PMID: 33052348; 33052332). Alternative splicing results in three transcript variants encoding distinct isoforms.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84510-4-RR (TMEM41B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.