

For Research Use Only

DCIR/CLEC4A Recombinant antibody, PBS Only (Capture)

Catalog Number: 84515-4-PBS



Basic Information

Catalog Number: 84515-4-PBS	GenBank Accession Number: NM_016184.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 50856	CloneNo.: 241854G1
Source: Rabbit	UNIPROT ID: Q9UMR7-1	
Isotype: IgG	Full Name: C-type lectin domain family 4, member A	
	Calculated MW: 28kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01361-1, DCIR/CLEC4A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84515-4-PBS. Detection antibody: 84515-1-PBS. Standard: Eg1350. Range: 0.469-60 ng/mL.