

For Research Use Only

LBR Recombinant antibody

Catalog Number: 84522-7-RR



Basic Information

Catalog Number: 84522-7-RR	GenBank Accession Number: BC020079	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 3930	CloneNo.: 241864G10
Source: Rabbit	UNIPROT ID: Q14739	Recommended Dilutions: IF/ICC 1:200-1:800
Isotype: IgG	Full Name: lamin B receptor	
Immunogen Catalog Number: AG3088	Calculated MW: 615 aa, 71 kDa	

Applications

Tested Applications: IF/ICC, ELISA	Positive Controls: IF/ICC : HeLa cells, HepG2 cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using LBR antibody (84522-7-RR, Clone: 241864G10) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using LBR antibody (84522-7-RR, Clone: 241864G10) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).