

For Research Use Only

Rat CD4 Recombinant antibody, PBS Only (Capture)

Catalog Number:84530-3-PBS



Basic Information

Catalog Number:

84530-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_012705.1

GeneID (NCBI):

24932

UNIPROT ID:

P05540

Full Name:

Cd4 molecule

Calculated MW:

51kDa

Purification Method:

Protein A purification

CloneNo.:

241955F10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01387-1, RAT CD4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84530-3-PBS. Detection antibody: 84530-2-PBS. Standard: Eg2025. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01387-2, RAT CD4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84530-3-PBS. Detection antibody: 84530-1-PBS. Standard: Eg2025. Range: 0.313-40 ng/mL