

For Research Use Only

IL1RN Recombinant antibody, PBS Only (Capture)

Catalog Number: 84532-1-PBS



Basic Information

Catalog Number: 84532-1-PBS	GenBank Accession Number: BC009745	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3557	CloneNo.: 241879D3
Source: Rabbit	UNIPROT ID: P18510	
Isotype: IgG	Full Name: interleukin 1 receptor antagonist	
Immunogen Catalog Number: AG1277	Calculated MW: 20 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01389-2, IL-1ra/IL-1F3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84532-1-PBS. Detection antibody: 84532-3-PBS. Standard: Ag1277. Range: 0.625-80 ng/mL.