

For Research Use Only

MDM4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84534-1-PBS



Basic Information

Catalog Number: 84534-1-PBS	GenBank Accession Number: BC067299	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4194	CloneNo.: 241936C10
Source: Rabbit	UNIPROT ID: O15151	
Isotype: IgG	Full Name: Mdm4 p53 binding protein homolog (mouse)	
Immunogen Catalog Number: AG26678	Calculated MW: 490 aa, 55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01394-2, MDMX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84534-1-PBS. Detection antibody: 84534-2-PBS. Standard: Ag26678. Range: 0.156-20 ng/mL