

For Research Use Only

MDMX/MDM4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84534-7-PBS



Basic Information

Catalog Number: 84534-7-PBS	GenBank Accession Number: BC067299	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4194	CloneNo.: 241936G5
Source: Rabbit	UNIPROT ID: O15151	
Isotype: IgG	Full Name: Mdm4 p53 binding protein homolog (mouse)	
Immunogen Catalog Number: AG26678	Calculated MW: 490 aa, 55 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01394-3, Human MDMX/MDM4 Recombinant Matched Antibody Pair - PBS only. 84534-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag26678. 84534-4-PBS was HRP conjugated as the detection antibody. Range: 19.5-1250 pg/mL



The mean MDMX/MDM4 concentration was determined to be 6,549.3 pg/mL in MCF-10A cell extract based on a 1.3 mg/mL extract load.