

For Research Use Only

FGF19 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84536-1-PBS



Basic Information

Catalog Number: 84536-1-PBS	GenBank Accession Number: BC017664	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9965	CloneNo.: 241740A12
Source: Rabbit	UNIPROT ID: O95750	
Isotype: IgG	Full Name: fibroblast growth factor 19	
Immunogen Catalog Number: HZ-1330	Calculated MW: 216 aa, 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01398-2, FGF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84536-1-PBS. Detection antibody: 84536-4-PBS. Standard: HZ-1330. Range: 0.078-20 ng/mL