

For Research Use Only

CD43 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84539-3-PBS



Basic Information

Catalog Number: 84539-3-PBS	GenBank Accession Number: NM_003123.6	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6693	CloneNo.: 242039B7
Source: Rabbit	ENSEMBL Gene ID: ENSG00000197471	
Isotype: IgG	UNIPROT ID: P16150	
	Full Name: sialophorin	
	Calculated MW: 40kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01360-2, CD43 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84539-4-PBS. Detection antibody: 84539-3-PBS. Standard: Eg1761. Range: 0.781-100 ng/mL.