

For Research Use Only

CCR8 Recombinant antibody, PBS Only

Catalog Number: 84545-1-PBS



Basic Information

Catalog Number: 84545-1-PBS	GenBank Accession Number: BC074743	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1237	CloneNo.: 241663E4
Source: Rabbit	UNIPROT ID: P51685	
Isotype: IgG	Full Name: chemokine (C-C motif) receptor 8	
	Calculated MW: 355 aa, 41 kDa	
	Observed MW: 41 kDa	

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 84545-1-RR (CCR8 antibody) at dilution of 1:800 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 84545-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 84545-1-RR (CCR8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84545-1-PBS in a different storage buffer formulation.