

For Research Use Only

DOCK2-Specific Recombinant antibody

Catalog Number: 84548-1-RR



Basic Information

Catalog Number:

84548-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004946

GeneID (NCBI):

1794

UNIPROT ID:

Q92608

Full Name:

dedicator of cytokinesis 2

Calculated MW:

212 kDa

Purification Method:

Protein A purification

CloneNo.:

241350D8

Recommended Dilutions:

IF/ICC 1:125-1:500

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : Jurkat cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using DOCK2-Specific antibody (84548-1-RR, Clone: 241350D8) at dilution of 1:250 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).