

For Research Use Only

Mouse CD21/CD35 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84565-1-PBS



Basic Information

Catalog Number: 84565-1-PBS	GenBank Accession Number: AAA37448.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 12902	CloneNo.: 241957E1
Source: Rabbit	UNIPROT ID: P19070	
Isotype: IgG	Full Name: complement receptor 2	
	Calculated MW: 113kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01407-2, MOUSE CD21 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84565-3-PBS. Detection antibody: 84565-1-PBS. Standard: Eg1396. Range: 0.781-100 ng/mL.