

For Research Use Only

LILRB3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84566-2-PBS



Basic Information

Catalog Number:

84566-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC112198

GeneID (NCBI):

11025

UNIPROT ID:

O75022

Full Name:

leukocyte immunoglobulin-like
receptor, subfamily B (with TM and
ITIM domains), member 3

Calculated MW:

631 aa, 69 kDa

Purification Method:

Protein A purification

CloneNo.:

241894E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01408-1, LILRB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84566-1-PBS. Detection antibody: 84566-2-PBS. Standard: Eg2084. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01408-2, LILRB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84566-3-PBS. Detection antibody: 84566-2-PBS. Standard: Eg2084. Range: 0.625-80 ng/mL