

For Research Use Only

PIGK Recombinant antibody, PBS Only (Capture)

Catalog Number: 84568-2-PBS



Basic Information

Catalog Number: 84568-2-PBS	GenBank Accession Number: BC026186	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10026	CloneNo.: 241989E6
Source: Rabbit	UNIPROT ID: Q92643	
Isotype: IgG	Full Name: phosphatidylinositol glycan anchor biosynthesis, class K	
Immunogen Catalog Number: AG3497	Calculated MW: 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01411-1, PIGK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84568-2-PBS. Detection antibody: 84568-1-PBS. Standard: Ag3497. Range: 0.781-100 ng/mL.