

For Research Use Only

Cortactin Recombinant antibody, PBS Only (Capture)

Catalog Number:84569-4-PBS



Basic Information

Catalog Number: 84569-4-PBS	GenBank Accession Number: BC008799	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2017	CloneNo.: 241987H3
Source: Rabbit	UNIPROT ID: Q14247	
Isotype: IgG	Full Name: cortactin	
Immunogen Catalog Number: AG1940	Calculated MW: 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01412-1, Cortactin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84569-4-PBS. Detection antibody: 84569-3-PBS. Standard: Ag1940. Range: 0.313-40 ng/mL.