

For Research Use Only

# RGNEF Recombinant antibody, PBS Only (Detector)

Catalog Number: 84570-1-PBS



## Basic Information

|   |   |   |
|---|---|---|
| <b>Catalog Number:</b><br>84570-1-PBS       | <b>GenBank Accession Number:</b><br>BC012946                | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>64283                              | <b>CloneNo.:</b><br>241975C10                         |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q8N1W1                                |   |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>Rho-guanine nucleotide exchange factor |   |
| <b>Immunogen Catalog Number:</b><br>AG10915 | <b>Calculated MW:</b><br>192 kDa                            |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01414-1, RGNLF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84570-2-PBS. Detection antibody: 84570-1-PBS. Standard: Ag10915. Range: 0.469-60 ng/mL