

For Research Use Only

WNT3A Recombinant antibody, PBS Only (Detector)

Catalog Number: 84577-1-PBS



Basic Information

Catalog Number:

84577-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25051

GenBank Accession Number:

BC103921

GeneID (NCBI):

89780

UNIPROT ID:

P56704

Full Name:

wingless-type MMTV integration site family, member 3A

Calculated MW:

352 aa, 39 kDa

Purification Method:

Protein A purification

CloneNo.:

241993A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01423-1, WNT3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84577-3-PBS. Detection antibody: 84577-1-PBS. Standard: Ag25051. Range: 0.781-100 ng/mL.