

For Research Use Only

GZMA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84583-2-PBS



Basic Information

Catalog Number: 84583-2-PBS	GenBank Accession Number: BC015739	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3001	CloneNo.: 241979E1
Source: Rabbit	UNIPROT ID: P12544	
Isotype: IgG	Full Name: granzyme A (granzyme 1, cytotoxic T-lymphocyte-associated serine esterase 3)	
Immunogen Catalog Number: AG1818	Calculated MW: 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01417-1, Granzyme A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84583-2-PBS. Detection antibody: 84583-1-PBS. Standard: Ag1818. Range: 0.625-80 ng/mL.