

For Research Use Only

BCAP31 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84584-3-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 84584-3-PBS | GenBank Accession Number: BC014323 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 10134 | CloneNo.: 241978F9 |
| Source: Rabbit | UNIPROT ID: P51572 | |
| Isotype: IgG | Full Name: B-cell receptor-associated protein 31 | |
| Immunogen Catalog Number: AG1668 | Calculated MW: 28 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01418-2, BCAP31 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84584-3-PBS. Detection antibody: 84584-4-PBS. Standard: Ag1668. Range: 0.625-80 ng/mL.