

For Research Use Only

DCX Recombinant antibody, PBS Only (Detector)

Catalog Number: 84595-2-PBS



Basic Information

Catalog Number:	84595-2-PBS	GenBank Accession Number:	BC027925	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	1641	CloneNo.:	241985E7
Source:	Rabbit	UNIPROT ID:	O43602		
Isotype:	IgG	Full Name:	doublecortin		
Immunogen Catalog Number:	AG4945	Calculated MW:	402 aa, 45 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01427-1, DCX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84595-3-PBS. Detection antibody: 84595-2-PBS. Standard: Ag4945. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01427-2, DCX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84595-1-PBS. Detection antibody: 84595-2-PBS. Standard: Ag4945. Range: 0.781-100 ng/mL.