

For Research Use Only

CD46 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84599-3-PBS



Basic Information

Catalog Number:

84599-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_172351.3

GeneID (NCBI):

4179

UNIPROT ID:

P15529-11

Full Name:

CD46 molecule, complement
regulatory protein

Calculated MW:

42kDa

Purification Method:

Protein A purification

CloneNo.:

241764E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01445-1, CD46 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84599-3-PBS. Detection antibody: 84599-2-PBS. Standard: Eg1770. Range: 0.469-60 ng/mL