

For Research Use Only

PZP Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84605-1-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 84605-1-PBS | GenBank Accession Number: BC111756 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 5858 | CloneNo.: 242114C11 |
| Source: Rabbit | Full Name: pregnancy-zone protein | |
| Isotype: IgG | Calculated MW: 1482 aa, 164 kDa | |
| Immunogen Catalog Number: AG16471 | | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01425-1, Pregnancy Zone Protein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84605-1-PBS. Detection antibody: 84605-3-PBS. Standard: Ag16471. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01425-2, Pregnancy Zone Protein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84605-2-PBS. Detection antibody: 84605-1-PBS. Standard: Ag16471. Range: 0.781-100 ng/mL