

For Research Use Only

PZP Recombinant antibody, PBS Only (Detector)

Catalog Number: 84605-3-PBS



Basic Information

Catalog Number: 84605-3-PBS	GenBank Accession Number: BC111756	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5858	CloneNo.: 242114G4
Source: Rabbit	Full Name: pregnancy-zone protein	
Isotype: IgG	Calculated MW: 1482 aa, 164 kDa	
Immunogen Catalog Number: AG16471		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01425-1, Pregnancy Zone Protein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84605-1-PBS. Detection antibody: 84605-3-PBS. Standard: Ag16471. Range: 0.781-100 ng/mL