

For Research Use Only

DPP4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84606-2-PBS



Basic Information

Catalog Number: 84606-2-PBS	GenBank Accession Number: BC013329	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1803	CloneNo.: 242045H4
Source: Rabbit	UNIPROT ID: P27487	
Isotype: IgG	Full Name: dipeptidyl-peptidase 4	
	Calculated MW: 88 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01439-1, DPP4/CD26 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84606-2-PBS. Detection antibody: 84606-3-PBS. Standard: Eg0157. Range: 0.781-50 ng/mL.