

For Research Use Only

PPM2C Recombinant antibody, PBS Only (Detector)

Catalog Number: 84612-2-PBS



Basic Information

Catalog Number: 84612-2-PBS	GenBank Accession Number: BC047619	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54704	CloneNo.: 242090D7
Source: Rabbit	UNIPROT ID: Q9P0J1	
Isotype: IgG	Full Name: protein phosphatase 2C, magnesium-dependent, catalytic subunit	
Immunogen Catalog Number: AG15453	Calculated MW: 537 aa, 61 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01450-1, PPM2C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84612-1-PBS. Detection antibody: 84612-2-PBS. Standard: Ag15453. Range: 1.25-80 ng/mL