

For Research Use Only

UTX Recombinant antibody, PBS Only (Capture)

Catalog Number: 84614-1-PBS



Basic Information

Catalog Number: 84614-1-PBS	GenBank Accession Number: BC093868	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7403	CloneNo.: 242033A9
Source: Rabbit	UNIPROT ID: O15550	
Isotype: IgG	Full Name: ubiquitously transcribed tetra-tricopeptide repeat, X chromosome	
Immunogen Catalog Number: AG36754	Calculated MW: 1401 aa, 154 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01453-2, UTX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84614-1-PBS. Detection antibody: 84614-2-PBS. Standard: Ag36754. Range: 0.781-100 ng/mL.