

For Research Use Only

SETD5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84622-2-PBS



Basic Information

Catalog Number: 84622-2-PBS	GenBank Accession Number: BC020956	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55209	CloneNo.: 241903B6
Source: Rabbit	UNIPROT ID: Q9C0A6	
Isotype: IgG	Full Name: SET domain containing 5	
Immunogen Catalog Number: AG28929		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01432-1, SETD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84622-2-PBS. Detection antibody: 84622-1-PBS. Standard: Ag28929. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP01432-2, SETD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84622-2-PBS. Detection antibody: 84622-3-PBS. Standard: Ag28929. Range: 0.469-60 ng/mL