

For Research Use Only

GMCL1 Recombinant antibody, PBS Only (Capture)

Catalog Number:84626-1-PBS



Basic Information

| | | |
|--|---|---|
| Catalog Number: 84626-1-PBS | GenBank Accession Number: BC007420 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 64395 | CloneNo.: 242028B9 |
| Source: Rabbit | UNIPROT ID: Q961K5 | |
| Isotype: IgG | Full Name: germ cell-less homolog 1 (Drosophila) | |
| Immunogen Catalog Number: AG7929 | Calculated MW: 59 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01436-1, GMCL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84626-1-PBS. Detection antibody: 84626-3-PBS. Standard: Ag7929. Range: 0.781-100 ng/mL.