

For Research Use Only

ARNT, HIF1B Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84650-2-PBS



Basic Information

Catalog Number: 84650-2-PBS	GenBank Accession Number: BC060838	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 405	CloneNo.: 242111E3
Source: Rabbit	UNIPROT ID: P27540	
Isotype: IgG	Full Name: aryl hydrocarbon receptor nuclear translocator	
Immunogen Catalog Number: AG5507	Calculated MW: 416 aa, 46 kDa, 87 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01470-1, ARNT/HIF-1 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84650-1-PBS. Detection antibody: 84650-2-PBS. Standard: Ag5507. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01470-2, ARNT/HIF-1 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84650-2-PBS. Detection antibody: 84650-3-PBS. Standard: Ag5507. Range: 0.781-100 ng/mL