

For Research Use Only

ATP8A2 Recombinant antibody, PBS Only (Capture)

Catalog Number:84653-1-PBS



Basic Information

Catalog Number:

84653-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36276

GenBank Accession Number:

NM_001313741

GeneID (NCBI):

51761

UNIPROT ID:

Q9NTI2

Full Name:

ATPase, aminophospholipid transporter-like, class I, type 8A, member 2

Calculated MW:

133kDa

Purification Method:

Protein A purification

CloneNo.:

242105D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01468-1, ATP8A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84653-1-PBS. Detection antibody: 84653-3-PBS. Standard: Ag36276. Range: 1.25-80 ng/mL.