

For Research Use Only

ATP8A2 Recombinant antibody, PBS Only (Capture)

Catalog Number:84653-2-PBS



Basic Information

Catalog Number: 84653-2-PBS	GenBank Accession Number: NM_001313741	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51761	CloneNo.: 242105G7
Source: Rabbit	UNIPROT ID: Q9NTI2	
Isotype: IgG	Full Name: ATPase, aminophospholipid transporter-like, class I, type 8A, member 2	
Immunogen Catalog Number: AG36276	Calculated MW: 133kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01468-2, ATP8A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84653-2-PBS. Detection antibody: 84653-3-PBS. Standard: Ag36276. Range: 1.25-80 ng/mL.