

For Research Use Only

# ATP8A2 Recombinant antibody, PBS Only (Detector)

Catalog Number:84653-3-PBS



## Basic Information

Catalog Number:

84653-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36276

GenBank Accession Number:

NM\_001313741

GeneID (NCBI):

51761

UNIPROT ID:

Q9NTI2

Full Name:

ATPase, aminophospholipid transporter-like, class I, type 8A, member 2

Calculated MW:

133kDa

Purification Method:

Protein A purification

CloneNo.:

242105H2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01468-2, ATP8A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84653-2-PBS. Detection antibody: 84653-3-PBS. Standard: Ag36276. Range: 1.25-80 ng/mL.



Cytometric bead array standard curve of MP01468-1, ATP8A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84653-1-PBS. Detection antibody: 84653-3-PBS. Standard: Ag36276. Range: 1.25-80 ng/mL.