

For Research Use Only

CXCR3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84655-3-PBS



Basic Information

Catalog Number:

84655-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC034403

GeneID (NCBI):

2833

UNIPROT ID:

P49682

Full Name:

chemokine (C-X-C motif) receptor 3

Calculated MW:

368 aa, 41 kDa

Purification Method:

Protein A purification

CloneNo.:

241790H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01488-1, CXCR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84655-3-PBS. Detection antibody: 84655-4-PBS. Standard: Eg2269. Range: 0.313-40 ng/mL