

For Research Use Only

# CXCR3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84655-4-PBS



## Basic Information

|                                       |   |   |
|---------------------------------------|---|---|
| <b>Catalog Number:</b><br>84655-4-PBS | <b>GenBank Accession Number:</b><br>BC034403            | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1 mg/ml               | <b>GeneID (NCBI):</b><br>2833                           | <b>CloneNo.:</b><br>241790H7                          |
| <b>Source:</b><br>Rabbit              | <b>UNIPROT ID:</b><br>P49682                            |   |
| <b>Isotype:</b><br>IgG                | <b>Full Name:</b><br>chemokine (C-X-C motif) receptor 3 |   |
|                                       | <b>Calculated MW:</b><br>368 aa, 41 kDa                 |   |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01488-1, CXCR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84655-3-PBS. Detection antibody: 84655-4-PBS. Standard: Eg2269. Range: 0.313-40 ng/mL.