

For Research Use Only

PSPC1 Recombinant antibody, PBS Only

Catalog Number: 84657-5-PBS



Basic Information

Catalog Number: 84657-5-PBS	GenBank Accession Number: BC014184	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55269	CloneNo.: 242158C1
Source: Rabbit	UNIPROT ID: Q8WXF1	
Isotype: IgG	Full Name: paraspeckle component 1	
Immunogen Catalog Number: AG10136	Calculated MW: 523 aa, 59 kDa	

Applications

Tested Applications:
IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using PSPC1 antibody (84657-5-RR, Clone: 242158C1) at dilution of 1:1000 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84657-5-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using PSPC1 antibody (84657-5-RR, Clone: 242158C1) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red). This data was developed using the same antibody clone with 84657-5-PBS in a different storage buffer formulation.