

For Research Use Only

# LAMP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84658-3-PBS



## Basic Information

Catalog Number:

84658-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005561.4

GeneID (NCBI):

3916

ENSEMBL Gene ID:

ENSG00000185896

UNIPROT ID:

P11279-1

Full Name:

lysosomal-associated membrane  
protein 1

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

242116E6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01489-2, LAMP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84658-3-PBS. Detection antibody: 84658-4-PBS. Standard: Eg2611. Range: 0.781-100 ng/mL.