

For Research Use Only

SIRP beta 1/CD172b Recombinant antibody, PBS Only (Detector)

Catalog Number: 84667-3-PBS



Basic Information

Catalog Number: 84667-3-PBS	GenBank Accession Number: NM_006065.5	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10326	CloneNo.: 242115D1
Source: Rabbit	UNIPROT ID: O00241-1	
Isotype: IgG	Full Name: signal-regulatory protein beta 1	
	Calculated MW: 43kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01478-2, SIRPB1/CD172b Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84667-1-PBS. Detection antibody: 84667-3-PBS. Standard: Eg1668. Range: 0.313-40 ng/mL.