

For Research Use Only

GANC Recombinant antibody, PBS Only (Detector)

Catalog Number: 84676-1-PBS



Basic Information

Catalog Number: 84676-1-PBS	GenBank Accession Number: BC059406	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2595	CloneNo.: 242060D2
Source: Rabbit	UNIPROT ID: Q8TET4	
Isotype: IgG	Full Name: glucosidase, alpha; neutral C	
Immunogen Catalog Number: AG35067	Calculated MW: 104 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01464-1, GANC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84676-2-PBS. Detection antibody: 84676-1-PBS. Standard:null. Range: null.