

For Research Use Only

# CDC40 Recombinant antibody, PBS Only (Detector)

Catalog Number:84679-1-PBS



## Basic Information

Catalog Number:

84679-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36151

GenBank Accession Number:

BC126114

GeneID (NCBI):

51362

UNIPROT ID:

O60508

Full Name:

cell division cycle 40 homolog (S.

cerevisiae)

Purification Method:

Protein A purification

CloneNo.:

242034G6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01501-2, CDC40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84679-3-PBS. Detection antibody: 84679-1-PBS. Standard: Ag36151. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01501-1, CDC40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84679-2-PBS. Detection antibody: 84679-1-PBS. Standard: Ag36151. Range: 0.313-40 ng/mL