

For Research Use Only

Jagged1 Recombinant antibody

Catalog Number: 84694-1-RR



Basic Information

Catalog Number:

84694-1-RR

Size:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC098393

GeneID (NCBI):

182

UNIPROT ID:

P78504

Full Name:

jagged 1 (Alagille syndrome)

Calculated MW:

1218 aa, 134 kDa

Observed MW:

150 kDa

Purification Method:

Protein A purification

CloneNo.:

241412F1

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells,

Background Information

Jagged1 is a type I, membrane- anchored, multidomain protein and a ligand of the Notch receptor. The receptor-ligand interaction can trigger the Notch signaling and activate transcription factors that play key roles in cell differentiation and morphogenesis. The predicted MW of Jagged1 is 134 kDa, while higher MW around 150-180 kDa had been observed in literatures (PMID: 9716576, 9596653).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 84694-1-RR (Jagged1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84694-1-RR against Human Jagged1 were performed. The affinity constant is 0.688 nM.