

For Research Use Only

# CYP51A1 Recombinant antibody, PBS Only (Capture)

Catalog Number:84705-2-PBS



## Basic Information

Catalog Number:

84705-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4235

GenBank Accession Number:

BC032322

GeneID (NCBI):

1595

UNIPROT ID:

Q16850

Full Name:

cytochrome P450, family 51, subfamily A, polypeptide 1

Calculated MW:

509 aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

242166D6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01526-2, CYP51A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84705-2-PBS. Detection antibody: 84705-3-PBS. Standard: Ag4235. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01526-1, CYP51A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84705-2-PBS. Detection antibody: 84705-1-PBS. Standard: Ag4235. Range: 0.781-100 ng/mL